

Practical Problems

Safety Regulations

The following regulations apply to all laboratories used for the Olympiad. Participating students must be well acquainted with these regulations and should study them seriously. These rules will be strictly followed in the 33rd IChO practical examination. Students who break any of these rules will be given only one warning before they are disqualified from the practical examination.

If any questions arise concerning safety procedures during the practical examination, students should not hesitate to ask the nearest instructor for directions.

All students are required to sign a statement agreeing that they have read and understood the rules prior to the examination.

Rules for personal safety

- a. For eye protection, safety goggles must be worn in the laboratories at all times. If the student wears contact lenses, full protection goggles, which provide total seal around eyes, must be worn. All students are requested to bring their safety goggles, but we shall have some in reserve.
- b. A long sleeved, knee length laboratory coat is recommended. Long pants and closed-toed shoes must be worn for individual safety. Loose clothing, open style shoes and sandals are prohibited. Long hair must be contained. Each student will have to get her/his own necessary items for herself/himself.
- c. Prior to the examination, the demonstrator-in-charge will check all protective equipments to ensure that they are in order.
- d. Pipetting by mouth is strictly forbidden.
- e. Eating, drinking or smoking in the laboratory is strictly prohibited.

Accidents and first aid

In any chemistry laboratory, accidents can take place due to spillage of chemicals, broken glasswares and fire. Any injury, illness, or incident, however minor, must be reported to the instructor immediately so that proper corrective action can be taken up.

a. Chemicals: Every chemical in the laboratory must be handled with utmost care. Chemicals can be corrosive, flammable or poisonous. Each student should read the safety notes related to the chemicals given in the task before handling them. The following general precautions must be always followed in the laboratory :

- ◆ Chemicals should never be tasted. Use pipette bulbs or pipette fillers all the time.
- ◆ Spillage on the skin: For any spillage of chemicals, the first step is to flush the skin under cold tap water for 10 to 15 minutes and then seek first aid/or medical help as appropriate. Organic materials tend to get absorbed on the skin, so wash the skin with warm water and soap, after cleaning it with cold water. Contaminated clothing should be removed at the earliest.
- ◆ Chemicals in the eye: The proper use of safety goggles will reduce the risk of any eye injury. Even so, if there is any splash of chemicals into the eyes, wash your eyes with cold water for 15 minutes and then look for appropriate medical attention.

b. Fire: Many chemicals are flammable, and hence no open flames are permitted when such chemicals are in use. You should get familiar with the location of the nearest fire extinguisher and fire blanket.

c. Glassware: Glass is a very hard but brittle material, and can break under stress or strain. Please handle the glasswares very carefully. If breakage occurs it is essential that any particles or splinters, specially from the wounds, are removed at the earliest. The injuries must be inspected by the demonstrator-in-charge.

Please report and clean up any breakage of the glassware. Necessary replacements can be obtained from the instructor.

- d. **Waste Materials:** Do not dispose of chemicals in the sink. Please follow all disposal rules provided in the task notes. Waste collection containers will be provided wherever necessary.
- e. **Care of Benches and Apparatus:** Each student is responsible for her/his section of the bench. Any spillage of chemicals or water must be wiped immediately. Concentrated acid spills must be first neutralized with sodium bicarbonate and then washed with plenty of water. Your working area must be kept clean at all times. Chemicals spilled on the ground must be washed and broken glassware must be swept off immediately. Mops, brooms, dust-pans etc will be available from the preparation room.

Some important information regarding the 33rd IChO practical examination

- Time duration for the practical examination would be four and a half hours instead of five hours.
- The examination may consist of three **independent** experimental tasks. The time duration for each task may vary from one to one and a half hour.
- The examination will be conducted in two batches. Students No.1 and 2 from each team will be part of the first batch; students No.3 and 4 will be part of the second batch.
- Students of both batches will get a new set of apparatus for the examination.
- The apparatus for the examination will include both plasticware and glassware.
- The examination will not involve use of microscale apparatus.


Blank titration

Dilute the 25 mL of 1 M NaOH solution in a 250 mL standard flask using freshly boiled distilled water. Pipette out 25 mL of the diluted NaOH solution and titrate it against the HCl solution using phenol red as indicator until the colour changes from red to yellow.

Titration of sample aliquot

Weigh accurately about 1.5 g of the crushed tablet sample and transfer it quantitatively in a 250 mL beaker. Add 25 mL of 1 M NaOH solution with the help of pipette and swirl the content. Boil the mixture gently on a water bath for 15 min and then cool the solution. Transfer the solution to a 250 mL standard flask. Dilute the solution up to the mark with distilled water and mix well. Titrate 25 mL of the diluted solution against the standardised HCl solution using phenol red indicator until the colour changes from red to yellow.

- *Write down the appropriate chemical reaction for hydrolysis of acetyl salicylic acid.*
- *Calculate the percentage of aspirin in the sample.*

Problem 26 Synthesis of 1-phenyl-azo-2-naphthol (C₁₆H₁₂ON₂)**Reactions**

Coupling reaction


Weigh 0.75 g of powdered β -naphthol in a 50 mL beaker. Add 5 mL of 10% NaOH solution and 5 mL of distilled water to the beaker. Stir well with glass rod to obtain a clear solution. This solution is also cooled in an ice-bath to 0°C.

The ice cooled filtrate containing diazotised salt is added dropwise to the ice cooled solution of β -naphthol with constant stirring. At this stage, an orange-red dye will start precipitating. After the addition of the solution is complete, filter the dye using buchner funnel. Cold water is used for washing the precipitate. Dry the product and record the yield.

Determination of melting point

Recrystallise a small portion of the organic dye prepared using ethyl alcohol. Gently heat the solution in a water bath (careful!) to dissolve the dye. Filter the hot solution. Cool the filtrate and filter the recrystallised product using Buchner funnel and suction.

- *Record the weight of the crude product*
- *Record the melting point of the recrystallised product.*

Problem 27 Determination of calcium in a sample solution**Reaction****Chemical and solutions**

- Sample solution containing calcium **R** : 36 **S** : 22, 24
(prepared from A.R. grade CaCl_2)
- Patton and Reeders indicator

- KOH solution. **R** : 35 **S** : 26, 37, 39, 45
- EDTA disodium salt **R** : 36, 37, 38 **S** : 26, 36

Preparation of 0.01 M EDTA:

Weigh 1.861 g of AR grade Na₂EDTA and quantitatively transfer the same to 500 mL volumetric flask. Add distilled water to the flask to dissolve Na₂EDTA and make up the solution to 500 mL mark with distilled water.

Procedure

Dilute the given sample solution to 100 mL in a 100 mL volumetric flask using distilled water. Pipette out 25 mL of the diluted sample solution in a clean conical flask. Add 25 mL of distilled water and adjust the pH using freshly prepared KOH solution to 12. Check the pH with pH paper. Add a pinch of solid indicator and titrate with Na₂EDTA solution till the colour changes from wine red to blue.

➤ **Calculate the amount of calcium in mmoles in 100 mL of the diluted sample solution**

Problem 28 Estimation of methyl ketone by back titration

Methyl ketones like acetone can be estimated by iodinating with excess of standard iodine in an alkaline medium. The unreacted iodine is then back titrated with standard sodium thiosulphate solution.

Chemicals and solutions

- 0.1N Iodine solution **R** : 20, 21 **S** : 23, 25
- 0.1N NaOH **R** : 35 **S** : 2, 26, 37, 39
- Concentrated HCl **R** : 34, 37 **S** : 26, 45
- 1 N H₂SO₄. **R** : 35 **S** : 2, 26, 30
- 0.1 M Na₂S₂O₃ **S** : 22, 24, 25

Preparation of 0.1 M Na₂S₂O₃:

Weigh 25 g of AR grade Na₂S₂O₃ and quantitatively transfer it to a 1 L volumetric flask. Prepare the solution using freshly boiled distilled water. Add 3 drops of chloroform while preparing the solution. Avoid exposure to light.

Preparation of 0.1 N I₂ solution

Dissolve 20 g of iodate-free potassium iodide in 30 - 40 mL of distilled water in a 1 L volumetric flask. Weigh 12.7 g iodine and quantitatively transfer to the concentrated potassium iodide solution. Shake the flask well until all the iodine dissolves and then dilute up to the mark with distilled water.

Procedure**Standardisation of Na₂S₂O₃**

Weigh out accurately 0.14 to 0.15 g of dry potassium iodate. Dissolve it in 25 mL of distilled and freshly boiled water and add 2 g of iodate free potassium iodide. Add 5 mL of 1N sulphuric acid. Titrate the liberated iodine with thiosulphate solution with constant shaking. When the colour of the solution is pale yellow add 200 mL of distilled water and 2 mL of starch indicator. Continue the titration until the colour changes from blue to colourless.

Determination of ketone

Weigh accurately 0.2 g of the given acetone sample in a clean 50 mL beaker and add minimum amount of distilled water. Transfer the acetone solution to a 250 mL standard volumetric flask. Add distilled water to the flask to prepare acetone solution in water and make up the solution to 250 mL mark with distilled water. Pipette out 10 mL of the acetone solution in a clean conical flask. Add 10 mL of 10% aqueous sodium hydroxide, and stopper the flask. Shake the flask for 10 min. At the end of 10 minutes, add 35 mL of 0.1 N Iodine solution from the burette. Swirl the content, preferably using magnetic stirrer for 5 minutes, and keep it standing for 15 minutes.

Yellow crystals of iodoform will appear. Acidify the solution with H₂SO₄ (check the pH with pH paper).

Preparation of 0.1 M Na₂S₂O₃

Weigh 25 g of AR grade Na₂S₂O₃ in a small beaker. Quantitatively transfer it to a 1 L volumetric flask. Prepare the solution using freshly boiled distilled water. Add 3 drops of chloroform while preparing the solution. Avoid exposure to light.

Standardisation of Na₂S₂O₃

Weigh out accurately 0.14 to 0.15 g of dry potassium iodate. Dissolve it in 25 mL of fresh, boiled distilled water and add 2 g of iodate free potassium iodide. Add 5 mL of 1N sulphuric acid. Titrate the liberated iodine with thiosulphate solution with constant shaking. When the colour of the solution is pale yellow add 200 mL of distilled water and 2 mL of starch indicator. Continue the titration until the colour changes from blue to colourless.

Procedure

Dissolve the given sample of phenol to 250 mL with distilled water. Take 25 mL of the phenol solution into 250 mL stoppered conical flask. Add 25 mL of standard potassium bromate solution and 0.5 g of potassium bromide. Add 5 mL of 3M sulphuric acid. Stopper the flask immediately. Mix the reagents and let them stand for 15 min (avoid exposure to light). Then, add 2.5 g of potassium iodide rapidly. Re-stopper the flask immediately and swirl the contents of the flask to dissolve the solid.

Titrate the liberated iodine with standard 0.1M Na₂S₂O₃ from the burette using starch indicator.

➤ ***Calculate the amount of phenol per 250 mL of the solution.***

Problem 30 Determination of amount of Fe (III) present in the given sample

Fe (III) in the sample solution is first reduced to Fe (II) in HCl medium using stannous chloride. Excess of stannous chloride is oxidized by addition of mercury (II) chloride. The Fe(II) is then titrated with standard potassium dichromate solution.

Chemicals and solutions

- Sample solution **R** : 36, 38 **S** : 26, 36
- 0.1N K₂Cr₂O₇ solution **R** : 45, 36, 37, 38, 43 **S** : 53, 22, 28
- Equimolar H₂SO₄ & **R** : 35 **S** : 2, 26, 30
 H₃PO₄ acid mixture **R** : 34 **S** : 26, 45
- Conc. HCl **R** : 34, 37 **S** : 26, 45
- 5% HgCl₂ **R** : 26, 27, 28 **S** : 13, 28, 45
- 3% SnCl₂ solutions **R**: 22, 36, 37, 38 **S** : 26, 36,
- Diphenylamine indicator. **R** : 23, 24, 25, 33 **S** : 28, 36, 37, 45

Note : NH₄Fe(SO₄)₂.12H₂O is used to prepare the sample solution

Preparation of 0.1N K₂Cr₂O₇ solution

Weigh accurately 1.225 g of pure K₂Cr₂O₇ and transfer it to a 250 mL volumetric flask. Prepare the solution using distilled water.

Procedure:

Dilute the given Fe(III) sample solution to 100 mL using the standard volumetric flask. Take 10 mL of the diluted sample solution in a clean conical flask. Add 2 mL of concentrated HCl and boil the solution. To the hot solution, add SnCl₂ solution dropwise till the reaction mixture becomes colourless. Add 2 - 3 drops of SnCl₂ in excess.

Cool the solution under tap water. Add 2 to 3 mL of HgCl₂ solution at once. A white precipitate is obtained at this stage. (If grey precipitate is obtained, reject the sample and start again.)

Add 2 to 3 mL of the acid mixture and 1 drop of the diphenylamine indicator and titrate it against K₂Cr₂O₇ solution. Continue the titration until a colour change from colourless to permanent blue or violet is observed.

- *Write down the appropriate chemical reactions .*
- *Calculate the amounts of Fe (III) and $\text{NH}_4\text{Fe}(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$ per 100 mL of the sample solution.*

